

Product Datasheet

TP53BP1 (Phospho-T543) antibody (orb544417)

Description	Rabbit polyclonal antibody to TP53BP1 (Phospho-T543)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T543 of human TP53BP1 (NP_001135451.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q12888
Entrez	7158
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.