

Product Datasheet

PKM (Phospho-Y105) antibody (orb544426)

Description	Rabbit polyclonal antibody to PKM (Phospho-Y105)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y105 of human PKM
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P14618
Entrez	5315
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.