

Product Datasheet

GYS1 (Phospho-S641) antibody (orb544428)

Description	Rabbit polyclonal antibody to GYS1 (Phospho-S641)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S641 of human Glycogen synthase 1 (GYS1)(NP_002094.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	100kDa
Uniprot ID	P13807
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.