


Product Datasheet

MAPT (Phospho-T498) antibody (orb544438)

Description	Rabbit polyclonal antibody to MAPT (Phospho-T498)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic phosphorylated peptide around T498 of human Tau(NP_058519.3?).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	50kDa
Uniprot ID	P10636
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeL...