

Product Datasheet

Canine Amyloid A Serum (SAA) protein [Out of stock] (orb544463)

Description	Recombinant Canine Amyloid A Serum (SAA) protein
Reactivity	Canine
Tested Applications	ELISA
Target	Amyloid A Serum (SAA) Canine Recombinant
Preservatives	No Preservative
Form/Appearance	Lyophilized
Note	For research use only
Application notes	Suitable for use in ELISA. May be used as a Control. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Source	E. coli
Expiration Date	6 months from date of receipt.