

Product Datasheet

Feline Amyloid A Serum (SAA) protein [Out of stock] (orb544465)

Description	Recombinant Feline Amyloid A Serum (SAA) protein
Species/Host	E. coli
Reactivity	Feline
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Amyloid A Serum (SAA) Feline Recombinant
Preservatives	No Preservative
Form/Appearance	Lyophilized
Note	For research use only
Application notes	Suitable for use in ELISA. May be used as a Control. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded
Clonality	Monoclonal
Expiration Date	6 months from date of receipt.