



Product Datasheet

GSK-3 Beta (phospho-Ser9) antibody (orb544567)


Description	GSK-3 Beta (phospho-Ser9) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GSK-3 Beta around the phosphorylation site of Ser9 (TT(p-S)FA)
Target	phospho-GSK-3 Beta (Ser9)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	47 kDa
Uniprot ID	P49841
Dilution Range	WB=1:500-2000, IHC-P=1:20-200
Expiration Date	12 months from date of receipt.



Flow cytometry of MCF7.Primary Antibody ...



IHC-P analysis of human breast using pho...



Western blot:MDA-MB-231(Human) Cell;293T...