



Product Datasheet

Peroxiredoxin 1 antibody (orb544576)

Description	Peroxiredoxin 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human Peroxiredoxin 1 protein, full length
Target	Peroxiredoxin 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	22 kDa
Uniprot ID	Q06830
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse liver using Peroxiredoxin 1 antibody



Western blot: HepG2 (Human) Cell; A549 (Human) Cell; L1210 (Mouse) Cell; Cerebrum (Mouse) Tissue