

Product Datasheet

MMP9 antibody (orb544577)

Description	MMP9 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human MMP9
Target	MMP9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	78 kDa
Uniprot ID	P14780
Expiration Date	12 months from date of receipt.