

Product Datasheet

Androgen Receptor(AR) antibody (orb544584)

Description	Androgen Receptor(AR) antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human human Androgen Receptor
Target	Androgen Receptor(AR)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	101 kDa
CAS Number	b3
Uniprot ID	P10275
Dilution Range	IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.