


Product Datasheet

BRCA1 (phospho-Ser1466) antibody (orb544693)

Description	BRCA1 (phospho-Ser1466) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human BRCA1 around the phosphorylation site of Ser1466 (PI(p-S)QN)
Target	Phospho-BRCA1 (Ser1466)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	208 kDa
Uniprot ID	P38398
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human
ovarian
cancer u...