



Product Datasheet

ANGPTL8 antibody (orb544742)

Description	ANGPTL8 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant mouse ANGPTL8 (16-198/198 aa)
Target	ANGPTL8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	20 kDa
Uniprot ID	Q6UXH0
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse liver
using ANGPTL8



Western
blot: Lane
1: Human
HepG2
cell; Lane...