

Product Datasheet

GSK-3 Beta antibody (HRP) (orb544772)

Description	GSK-3 Beta antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GSK-3 Beta
Target	GSK-3 Beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P49841
Dilution Range	WB=1:500-1000,IHC-P=1:100-500,IHC-F=1:200-500,ELISA=1:1000-5000
Expiration Date	12 months from date of receipt.