

Product Datasheet

Human SAE1 protein (orb54551)

Description	Recombinant human SAE1 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	N-terminal GST-tagged: N-terminal GST-tagged1-261AA: 1-346AFull Length : Full Length
Protein Sequence	MVEKEEAGGGISEEEAAQYDRQIRLWGLEAQKRLRASRVLLVGLKGLGAEIAKN LILAGVKGLTMLDHEQVTPEDPGAQFLIRTGSVGRNRAEASLERAQNLNPMVDV KVDTEIDIEKKPESFFTQFDAVCLTCCSRDVIVKVDQICHKNSIKFFTGDVFGYHGY TFANLGEHEFVEEKTAKVVSQGVEDGPDTKRAKLDSETTMVKKKVVFCPVKE ALEVDWSSEKAKAALKRTTSDYFLLQVLLKFRTDKGRDPSSDYEEDSELLQIR NDVLDLSLGISPDLLPEDFVRYCFSEMAPVCAVVGILAQEIVKALSQRDPPHNNF FFFDGMKNGNGIVECLGPK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	65.4 kDa
Uniprot ID	Q9UBE0
Expiration Date	6 months from date of receipt.