


Product Datasheet

ATP Synthase Microplate Assay Kit (orb545613)

Description	Detection and Quantification of ATP Synthase Activity.
Reactivity	All
Range	0.1 µmol/ml - 10 µmol/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve