

Product Datasheet

Granule Bound Starch Synthase Microplate Assay Kit (orb545615)

Description	Detection and Quantification of Granule Bound Starch Synthase Activity.
Reactivity	All
Range	4 µmol/L - 400 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9030-10-8
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.