


Product Datasheet

Albumin (BCG) Microplate Assay Kit (orb545619)

Description	Detection and Quantification of Albumin Content.
Reactivity	All
Range	1 mg/ml - 100 mg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Urine, Serum, Plasma, Other biological fluids
CAS Number	70024-90-7; 9006-59-1; 9048-46-8
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve