

Product Datasheet

Thiol Microplate Assay Kit (orb545623)

Description Detection and Quantification of Thiol Content.

Reactivity All

Range 0.02 mmol/L - 2 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Colorimetric

Sample Types Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

Hazard Information Non-Toxic

Expiration Date 12 months from date of receipt.



Representative
Standard
Curve