


Product Datasheet

Urea Microplate Assay Kit (orb545624)

Description	Detection and Quantification of Urea Content.
Reactivity	All
Range	5 mg/L - 250 mg/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Urine, Serum, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	57-13-6
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve