

Product Datasheet

Isocitrate Dehydrogenase Microplate Assay Kit (orb545631)

Description Detection and Quantification of Isocitrate Dehydrogenase Activity.

Reactivity All

Range 0.2 mmol/L - 20 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.

Note For research use only


Application notes Colorimetric

Sample Types Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 9028-48-2

Hazard Information Non-Toxic

Expiration Date 12 months from date of receipt.



Representative
Standard
Curve