



Product Datasheet

Urease Microplate Assay Kit (orb545636)

Description	Detection and Quantification of Urease Activity.
Reactivity	All
Range	0.05 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Urine, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9002-13-5
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...