


Product Datasheet

Ammonia/Ammonium Microplate Assay Kit (orb545637)

Description	Detection and Quantification of Ammonia/Ammonium Content.
Reactivity	All
Range	1 µmol/L - 1000 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Serum, Plasma, Urine, Saliva, Cell culture media, Tissue extracts, Cell lysate, Other biological fluids
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve