


Product Datasheet

Soluble Sugar Microplate Assay Kit (orb545638)

Description	Detection and Quantification of Soluble Sugar Content.
Reactivity	All
Range	2 µg/ml - 200 µg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Serum, Plasma, Cell culture media, Tissue extracts, Cell lysate, Other biological fluids
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve