


Product Datasheet

endo-beta-Mannanase Microplate Assay Kit (orb545639)

Description	Detection and Quantification of endo-beta-Mannanase Activity.
Reactivity	All
Range	1 mmol/L - 10 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Tissue extracts, Cell lysate, Other biological fluids
CAS Number	37288-54-3
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve