


Product Datasheet

beta-N-Acetylglucosaminidase Microplate Assay Kit (orb545640)

Description	Detection and Quantification of beta-N-Acetylglucosaminidase Activity.
Reactivity	All
Range	5 µmol/L - 500 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Serum, Plasma, Cell culture media, Urine, Other biological fluids
CAS Number	9012-33-3
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve