



Product Datasheet

CD164 antibody (orb545671)

Description	Mouse monoclonal antibody to CD164
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	Breast tumor cell line T-47D
Target	CD164
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml. Western blotting: Recommended dilution: 1-2 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q04900
Entrez	8763
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml. Western blotting: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...



Western
blotting
analysis of
human
CD164...