


Product Datasheet

CD167a antibody (orb545673)

Description	Mouse monoclonal antibody to CD167a
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	CD167a-transfected NIH-3T3 cells
Target	CD167a
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Uniprot ID	Q08345
Entrez	780
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...