



Product Datasheet

CD169_Sialoadhesin antibody (APC) (orb545675)


Description	Mouse monoclonal antibody against
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	human rhinovirus 14-infected monocyte-derived dendritic cells
Target	CD169
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	Q9BZZZ
Entrez	6614
Dilution Range	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation
of human
CD169
positive
CD11c...



Flow
cytometry
surface
staining
pattern ...