



Product Datasheet

CD263 antibody (APC) (orb545690)

Description	Mouse monoclonal antibody against CD263 conjugated to APC
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	TRAIL-R3-hlgGhc fusion protein
Target	CD263
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	O14798
Entrez	8794
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD263 positive cells...



Flow cytometry surface staining pattern ...