


Product Datasheet

CD305 antibody (APC) (orb545697)


Description

Mouse monoclonal antibody against CD305 conjugated


Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Activated NK cells and CD3- thymocytes
Target	CD305
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	Q6GTX8
Entrez	3903
Dilution Range	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD305
positive
CD19 ...



Flow
cytometry
surface
staining
pattern ...