



Product Datasheet

HLA-F antibody (orb545716)

Description	Mouse monoclonal antibody to HLA-F
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Inclusion body-derived HLA-F heavy chain
Target	HLA-F
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P30511
Entrez	3134
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human activated lymphocyte...



Flow cytometry surface staining pattern ...