



Product Datasheet

MICA/MICB antibody (orb545718)

Description	Mouse monoclonal antibody to MICA/MICB
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP
Immunogen	Transfected C1R cells expressing MICA
Target	MICA/MICB
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human HeLa cells (red-fill...)



Flow cytometry surface staining pattern ...