


Product Datasheet

Human GM-CSF R alpha Protein (orb545735)

Description	GM-CSF R alpha(Granulocyte-macrophage colony-stimulating factor receptor subunit alpha) is also know...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	GM-CSF R alpha
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 60.9 kDa. The protein migrates as 70-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human GM-CSF R alpha, Fc Tag (orb545735) is expressed from human 293 cells (HEK293). It contains AA Glu 23 - Gly 320 (Accession # P15509-1).
MW	60.9 kDa
NCBI	NP_006131.2
Expiration Date	6 months from date of receipt.



Loaded
Human GM-
CSF R alpha,
Fc Tag
(Cat...