

Product Datasheet

Human B7-H6 / NCR3LG1 Protein (orb545764)

Description	Human B7-H6 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Asp 25 - Ser...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	B7-H6
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.6 kDa. The protein migrates as 38-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP 001189368.1
Purity	90%
Source	Human B7-H6 Protein, His Tag (orb545764) is expressed from human 293 cells (HEK293). It contains AA Asp 25 - Ser 262 (Accession # Q68D85-1).
MW	28.6 kDa
Uniprot ID	Q68D85
NCBI	001189368.1
Expiration Date	6 months from date of receipt.