

Product Datasheet

Leupeptin hemisulfate (orb545887)

Description	Leupeptin hemisulfate
Target	Protease
Form/Appearance	powder
Storage	-20°C
Note	For research use only
Formula	$C_{20}H_{38}N_6O_4 \cdot 1/2H_2SO_4$
Purity	> 95%; peptide content > 85%
MW	475.6
CAS Number	103476-89-7
Solubility (25°C)	Soluble in Water (up to 50 mg/ml). Also soluble in ethanol and methanol
Expiration Date	12 months from date of receipt.