



Product Datasheet

THOC1 antibody (orb5460)

Description	THOC1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human HPR1/Nuclear Matrix Protein p84 (221-320/657 aa)
Target	Nuclear Matrix Protein p84
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72 kDa
Uniprot ID	Q96FV9
Expiration Date	12 months from date of receipt.



Western blot analysis of rat kidney lysa...



Paraffin-embedded human lung carcinoma f...