

Product Datasheet

GSTA1/GSTA2/GSTA3/GSTA4/GSTA5 Biotinylated Antibody (orb546312)

Description

GSTA1/GSTA2/GSTA3/GSTA4/GSTA5

Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	E. coli-derived humn GSTA1/A2/A3/A4/A5 recombinant protein (Position: A2-F222).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25722 MW
Uniprot ID	P08263
Dilution Range	ELISA Detection (Matched Antibody Pair), 0.1-0.5µg/ml
Expiration Date	12 months from date of receipt.