

Product Datasheet

ARID2 Antibody (orb546354)

Description
ARID2 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence of human ARID2 (DERRKGLAFLDELRFHHSR).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


200 kDa

Uniprot ID
Q68CP9
Dilution Range


Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml

Expiration Date


12 months from date of receipt.



WB analysis of ARID2 using anti-ARID2 an...



IHC analysis of ARID2 using anti-ARID2 a...



IHC analysis of ARID2 using anti-ARID2 a...