

Product Datasheet

RAB3GAP2 Antibody (orb546373)

Description

RAB3GAP2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

E.coli-derived human RAB3GAP2 recombinant protein (Position: M1-Q115).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

156 kDa

Uniprot ID


[Q9H2M9](#)

Dilution Range

Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Direct ELISA, 0.1-0.5µg/ml


Expiration Date

12 months from date of receipt.




Flow
Cytometry
analysis of
CACO-2 cells

...



WB analysis
of RAB3GAP2
using anti-
RAB3G...



WB analysis
of RAB3GAP2
using anti-
RAB3G...