

Product Datasheet

CPAMD8 Antibody (orb546399)

Description
CPAMD8 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

E.coli-derived human CPAMD8 recombinant protein
(Position: R58-D234).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of
500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After
reconstitution, at 4°C for one month. It can also be
aliquotted and stored frozen at -20°C for six months. Avoid
repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive
results. By Heat: Boiling the paraffin sections in 10mM
citrate buffer, pH6.0, for 20mins is required for the staining
of formalin/paraffin sections. Other applications have not
been tested. Optimal dilutions should be determined by
end users. Add 0.2ml of distilled water will yield a
concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW




207 kDa

Uniprot ID
Q8IZJ3
Dilution Range

Western blot, 0.1-0.5µg/ml Immunohistochemistry
(Paraffin-embedded Section), 0.5-1µg/ml Direct ELISA, 0.1-
0.5µg/ml

Expiration Date

12 months from date of receipt.


WB analysis
of CPAMD8
using anti-
CPAMD8 ...

IHC analysis
of CPAMD8
using anti-
CPAMD8...

IHC analysis
of CPAMD8
using anti-
CPAMD8...