

Product Datasheet

Transmembrane protein 199 TMEM199 Antibody (orb546400)

Description

Transmembrane protein 199 TMEM199 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

E.coli-derived human TMEM199 recombinant protein (Position: E20-H129).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

23 kDa

Uniprot ID

[Q8N511](#)


Dilution Range

Western blot, 0.1-0.5µg/ml Direct ELISA, 0.1-0.5µg/ml

Expiration Date

12 months from date of receipt.

1 2 3 4 5 6 7 8



WB
analysis of
TMEM199
using anti-
TMEM19...