


Product Datasheet

Pig porcine Activin A ELISA Kit (orb546404)

Description	Pig porcine Activin A ELISA Kit
Reactivity	Porcine
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: G311-S426
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva.
MW	47442 MW
Uniprot ID	A4D1W7
Sensitivity	12 pg/ml
Expiration Date	Please enquire.



Pig Activin A
PicoKine
ELISA Kit
standar...