


Product Datasheet

Monkey primate BDNF ELISA Kit (orb546410)

Description	Monkey primate BDNF ELISA Kit
Reactivity	Monkey
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H129-R247
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	27818 MW
Uniprot ID	P23560
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Monkey
Primate
BDNF
PicoKine
ELISA Kit s...