


Product Datasheet

Monkey primate BMP-5 ELISA Kit (orb546417)

Description	Monkey primate BMP-5 ELISA Kit
Reactivity	Monkey
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A317-H454
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and bone tissue.
MW	51737 MW
Uniprot ID	P22003
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Monkey
Primate
BMP-5
PicoKine
ELISA Kit ...