


Product Datasheet

Human Erythropoietin/EPO ELISA Kit (orb546427)

Description	Human Erythropoietin/EPO ELISA Kit
Reactivity	Human
Range	3.12 mIU/ml - 200 mIU/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: A28-R193
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P01588
Sensitivity	1mIU/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is around 60 mIU/m. For plasma, dilution ratio is 1:1 and concentration is 3 mIU/m.
Expiration Date	Please enquire.



Human EPO
PicoKine
ELISA Kit
Standard
Cu...