


Product Datasheet

Human FLT3/FLK2 ELISA Kit (orb546434)

Description	Human FLT3/FLK2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: N27-N541
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, lysates, tissue, serum and plasma (heparin, EDTA).
Uniprot ID	P36888
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 100 pg/ml.
Expiration Date	Please enquire.



Human
FLT3/FLK2
PicoKine
ELISA Kit
Stand...