


Product Datasheet

Pig porcine TGF Alpha ELISA Kit (orb546448)

Description	Pig porcine TGF Alpha ELISA Kit
Reactivity	Porcine
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V40-A89
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and Pig porcine milk.
MW	17006 MW
Uniprot ID	P01135
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Pig TGF
alpha
PicoKine
ELISA Kit
standar...