


Product Datasheet

Rabbit TGF Beta 1 ELISA Kit (orb546455)

Description	Rabbit TGF Beta 1 ELISA Kit
Reactivity	Rabbit
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: A279-S390
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (EDTA) and urine.
MW	44341 MW
Uniprot ID	P01137
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Rabbit TGF
beta 1
PicoKine
ELISA Kit
sta...