


Product Datasheet

Dog Canine Endothelin 1/EDN1 ELISA Kit (orb546475)

Description	Dog Canine Endothelin 1/EDN1 ELISA Kit
Reactivity	Canine
Range	3.9 pg/ml - 250 pg/ml
Tested Applications	ELISA
Immunogen	A peptide reflecting the target protein in human. Sequence: C53-W73
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	24425 MW
Uniprot ID	P12064, P13206, Q765Z4
Sensitivity	0.5 pg/ml
Expiration Date	Please enquire.



Dog
Endothelin
PicoKine
ELISA Kit
standa...