


Product Datasheet

Horse equine Endothelin 1/EDN1 ELISA Kit (orb546476)

Description	Horse equine Endothelin 1/EDN1 ELISA Kit
Reactivity	Equine
Range	3.9 pg/ml - 250 pg/ml
Tested Applications	ELISA
Immunogen	A peptide reflecting the target protein in horse equine. Sequence: C53-W73
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	24425 MW
Uniprot ID	F7APB1, F7AE96, Q867D0
Sensitivity	0.5 pg/ml
Expiration Date	Please enquire.



Horse
equine
Endothelin
PicoKine
ELISA K...